

# **MEDIUM-DUTY AMBULANCES**

**Crew-tailored Emergency Vehicles**



With a legacy and passion to deliver mission-critical durability, Wheeled Coach ambulances are born to perform under the most demanding conditions—and built to be as unbreakable as the spirit of those we serve. That's why they are truly ***Trusted by the Toughest.***

We understand every department has its own unique requirements when it comes to their emergency vehicles. That's why we build flexibility and superior performance into every ambulance we manufacture. Our Medium-Duty ambulances are built tough on a Ford, Chevy, Freightliner or International heavy-duty chassis to take on the most demanding conditions. Add more room, options or both to make these Medical Attack Vehicles, as they're known, exactly what you need for your life-saving missions.

## **TOP FEATURES**



**CORNER CAP LIGHTS**



**SCBA BRACKETS IN COMPARTMENT B**



**SAFETY-MINDED CABINET DESIGN**



**REV COMFORT CHILD SAFETY SEAT**



**COOL-BAR™**



**REAR DOOR GRAB BARS**

# SPECIFICATIONS

## Medium-Duty Model Specifications

	Chassis	Model	Length	Width	Headroom
<b>Ford</b>	F650	9170	170"	95"	72"
<b>Chevy</b>	MD4500/5500	9170	170"	95"	72"
<b>Freightliner</b>	Freightliner	9170	170"	95"	72"
<b>International</b>	International	9170	170"	95"	72"

## Standard Features

- |                                                                                                                                                                                                                                                                                                                                                                       |                                                                                                                                                                                                                                                                                                                                                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>• 20-year structural body warranty</li> <li>• Black Poly-urea coated thresholds</li> <li>• Running boards with Gator Grip surface</li> <li>• Full floating automotive style compartment door handles</li> <li>• All-aluminum modular body built in house</li> <li>• Lonplate II raised diamond plate vinyl flooring</li> </ul> | <ul style="list-style-type: none"> <li>• TrimMark OEM free-floating door handles</li> <li>• Large handicap style grab handles on side and rear doors with powder coated safety yellow coating</li> <li>• Aluminum console</li> <li>• Per4Max® seat belt system</li> <li>• Rear high-back captain's chair attendant seat with 4-point harness</li> <li>• Poly-urea coated compartments</li> </ul> | <ul style="list-style-type: none"> <li>• Reflective rub rails with multiple color choices</li> <li>• SafePASS™ emergency egress lock releases on side and rear patient compartment doors</li> <li>• Seamless upholstery</li> <li>• Single sheet seamless construction, no exposed door jambs</li> <li>• Stainless-steel curbside and rear door panels with reflective safety chevrons</li> </ul> |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## WHEELED COACH FEATURES

Wheeled Coach is an industry leader in product innovations that protect the occupants while providing long life and durability. Our Medium-Duty ambulances include these built-in features.

### Safeguard Box Frame™

Using advanced processes that ensure the highest quality, repeatable results and precision, we build a highly reinforced, but ultra-light roll cage superstructure with strategically positioned metal bars, creating one of the strongest, safest and most durable ambulance bodies around.

### AntiCor™ Protection System

Layer upon layer of anti-corrosion features that stand up to the most demanding conditions.

### ToughCoat™ Paint Process

Our unique nine-step process employs precision application standards centered around a revolutionary high-tech primer with adhesion and anti-corrosive properties never before achieved for a super tough, long-lasting finish.

### TrueForm™ Wall Construction

This precision assembly process incorporates our **StraightLine™ Laser-guided Fixture Alignment** to create accurate, square, repeatable frames that fit every time.

Our network of independent dealers, backed by the Wheeled Coach team and REV Finance, are ready as your first line of defense when you need support. Contact our sales team at **1-800-422-8206** or find your local dealer at: [wheeledcoach.com/ambulance-dealer-locator](http://wheeledcoach.com/ambulance-dealer-locator).

